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<http://www.iceotherm.com>

ICE-O-therm® PZ2

Congratulations on your purchase of the ICE-O-therm® PZ2 from Saringer Life Science Technologies Inc. 13231 Kennedy Road, Stouffville, ON, Canada, L4A 7X5. The patented ICE-O-therm® PZ2 technology offers temporary relief of minor aches and pains and muscle spasms by applying warm and cool stimuli to a specific area on the body. Please read these instructions carefully before using your device.

Getting the ICE-O-therm® ready to use

- The ICE-O-therm®PZ2 will run for 1 - 2 hours on a set of fully charged batteries.
- CAUTION:
 - Do not use 1 charged, and 1 partially charged battery in the unit at the same time. Always put 2 fully charged batteries into the unit.
- Unscrew the bottom of the ICE-O-therm®PZ2 unit (where the carrying strap is located).
- Insert 2 batteries (positive side towards the tip) into the unit, and screw the battery cap back on.
- Press the button on the bottom of the unit to turn it on (the LED lights will glow, and the white square pad will start to feel cool to the touch). The ICE-O-therm®PZ2 will take 4 - 5 minutes to reach its optimum operating temperature.
- Temperature range is between 5 and 50 degrees centigrade.

Using the ICE-O-therm®PZ2

- Once the ICE-O-therm®PZ2 has reached its optimum operating temperature (4 – 5 minutes), place the flat part of the ICE-O-therm®PZ2 head (with the white pad) on a small joint, or other area where you are experiencing pain.
- To achieve effective pain relief, try to maintain contact between the flat part of the ICE-O-therm® treatment head and the area you are experiencing pain.
- Some people find the sensation that the ICE-O-therm®PZ2 produces while in contact with skin to be quite intense. Rest assured that the ICE-O-therm® PZ2 cannot harm you if used correctly.
- Hold the unit in place and slowly move it around the affected area.
- A typical treatment will take 10 – 20 minutes, depending on the person and type of pain.
- Do not apply the unit for more than 20 minutes on 1 location.
- NOTE: When the batteries start to wear out, the white pad will not feel as cool to the touch. When this occurs, replace the batteries with fresh ones.

Precautions

- Do not use for more than 20 minutes on a given area at one time.
- Do not apply to broken or sensitive skin.
- If you suffer from hyperesthesia, please consult your physician before using the ICE-O-therm®PZ2.

Troubleshooting

- If you are having trouble with your device (LED lights do not come on, or you feel no warm / cold coming from the head), please try the following before contacting customer support:
 - Ensure that the unit has 2 fully charged batteries inside, and that they are inserted in the correct orientation.
 - Ensure that the battery cap is screwed on tightly.
 - If you are still experiencing trouble with your device, please contact us at: customerservice@saringer.com

Included in package

- ICE-O-therm®PZ2 unit with carrying strap

- ICE-O-therm®PZ2 carrying case
- 4 AA Alkaline batteries in carrying case

ICE-O-therm is not a sterile device. Caution should be exercised when sharing amongst multiple users. Wipe head with an alcohol swab.

WARNING: Use carefully. May cause burns. Do not use over sensitive skin areas or in the presence of poor circulation. The unattended use of this device by children or incapacitated persons may be dangerous.